

## SECONDARY ANTIBODIES COATED SURFACE: GOAT ANTI MOUSE IgG Fc $\gamma$ (Subclasses 1+2a+2b+3)

### TECHNICAL NOTE N. 37

#### Binding capacity and sensitivity test

1. Add 100  $\mu$ l of different concentrations of mouse IgG (from 0.025 to 4  $\mu$ g/ml) to the wells of goat anti mouse IgG coated plate and incubate for 60 minutes at room temperature
2. Empty the wells and wash with 0.1 M PBS pH 7.2, 0.05% Tween<sup>®</sup> 20 four times
3. Add 100  $\mu$ l /well of Goat anti-mouse IgG (H+ L)-HRP (Jackson ImmunoResearch code 115-035-003, diluted 1: 150.000) and incubate for 30 minutes at room temperature
4. Empty the wells and wash with 0.1 M PBS pH 7.2, 0.05% Tween<sup>®</sup> 20 four times
5. Add 100  $\mu$ l /well of TMB substrate solution and incubate 15 minutes at room temperature
6. Stop the substrate reaction by adding 100  $\mu$ l /well of sulphuric acid 0.3 N and read the optical density values at 450 nm

The data show that a plateau has got starting with an IgG mouse concentration of 1.0  $\mu$ g/ml.

This concentration means the well binding capacity we can express as:

- $\mu$ g/well = 0.1 (100 ng/well)
- pmol/well = 0.625 (this result is calculated considering the IgG M.W. = 160.000)

The microplate sensitivity was calculated as the lowest mouse IgG concentration higher than the mean optical density plus 5 S.D. of 0  $\mu$ g/ml mouse IgG concentration.

Our experiment gave the following results:

- 0  $\mu$ g/ml mouse IgG optical density mean (coming from 8 replicates) = 0.108
- standard deviation = 0.014
- mean + 5 S.D. = 0.178
- sensitivity = 0.012  $\mu$ g/well of mouse IgG

**Figure 1**

