

## PRODUCT DATA SHEET

### NO BINDING SURFACE

96 well plate

#### 1. Description

Code	Plate format	Colour	Volume per well
MC03F-NB	96 well solid plate	Clear	400 µl
MCW03F-NB	96 well solid plate	White	400 µl
MCB03F-NB	96 well solid plate	Black	400 µl

plate	Material	Polystyrene
	Well bottom	Flat
	Well shape	Round
Dimensions	SBS standard	

#### 2. Features

Surface	No Binding
Active Surface	300 µl
Uniformity	CV% < 10% with our standard test. For further details see Certificate of Quality
Packaging	Plates are packed in a transparent bag
Storage	Room Temperature. Once opened, place unused plate in a bag and store at room temperature.
Shelf life	2 years from date of manufacture
Sterile	No
Ready to use	Yes
Other information	All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit [www.biomat.it](http://www.biomat.it)