

## PRODUCT DATA SHEET

### HTS PLATES

96 well plate

#### 1. Description

Code	Plate format	Strips colour	Volume per well	Breakable
MC0F-96	96 well solid plate	Clear	400 µl	No
MCW0F-96	96 well solid plate	White	400 µl	No
MCB0F-96	96 well solid plate	Black	400 µl	No

Strip	Material	n/a
	Well bottom	Flat
	Well shape	Round
Frame	Colour	See table
	Material	Polystyrene
Dimensions	SBS standard	

#### 2. Features

Packaging	Qty/ pk: 25
	Qty/ cs: 200
Storage	Room temperature
Other information	All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

For detailed technical information visit [www.biomat.it](http://www.biomat.it)