

PRODUCT DATA SHEET

96 Deep Well Plate 1.2 ml

1. Description

Code	Plate format	Colour	Volume per well
DWRU1200	96 Deep Well plate	Clear	1.2 ml

Plate	Material	Polypropylene (high quality)
	Well shape	Round
	Well bottom	U bottom

2. Features

Packaging	Qty/ pk: 5
	Qty/ cs: 50
Storage	Room temperature in a dry place. Keep out of sunlight.
Other information	DNase and RNase free Good chemical resistance and heat-resistance.

Data Sheet subject to change without notice.

For detailed technical information visit www.biomat.it