

PRODUCT DATA SHEET

HIGH BINDING CAPACITY

96 well plate

1. Description

Code	Plate format	Strips colour	Volume per well
MC01F-HB8	96 well solid plate	Clear	400 µl
MCW01F-HB8	96 well solid plate	White	400 µl
MCB01F-HB8	96 well solid plate	Black	400 µl

plate	Material	Polystyrene
	Well bottom	Flat
	Well shape	Round
Dimensions	SBS standard	

2. Features

Surface	High binding
Binding capacity	~ 400 ng IgG/cm ²
Uniformity	CV% < 5% with our standard test. For further details see Certificate of Quality
Packaging	Qty/ pk: 25
	Qty/ cs: 200
Storage	Room temperature
Shelf life	5 years
Other information	All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it