

PRODUCT DATA SHEET

AMMINATED SURFACE

96 well plate

1. Description

Code	Plate format	Colour	Volume per well
MC02F-AM1	96 well solid plate	Clear	400 µl
MCW02F-AM1	96 well solid plate	White	400 µl
MCB02F-AM1	96 well solid plate	Black	400 µl

plate	Material	Polystyrene
	Well bottom	Flat
	Well shape	Round
Dimensions	SBS standard	

2. Features

Surface	Primary Amino Groups
Active Surface	200 µl
Packaging	Plates are packed in a transparent bag
Storage	Room Temperature. Once opened, place unused plates in a bag and store at room temperature.
Shelf life	2 years from date of manufacture
Sterile	No
Ready to use	Yes
Other information	All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it