

PRODUCT DATA SHEET

BARRIER BAGS

1. Description

Code	Description	Colour	Material	Dimension
AB-rs	Barrier bags for microplates, re-sealable with zip	Silver	see table	Internal: mm 115 x 175 External: mm 135 x 210

2. Features

Composition	12 µm Polyester (PET)
	8 µm Aluminium foil
	100 µm Linear low-density polyethylene (LLDPE)
Total thickness	102 µm (± 5%)
Mass	120 g/m ²
Water vapour transmission rate	0.01 g/m ² /24 h (38° - 95% RH)
Classification	MIL-B-131-H-TYPE 1-CLASS 1
Packaging	Q.ty/case: 1500
Other information	Re-sealable with zip. Zip is 25 mm from sealed end ± 0.3 mm

Data Sheet subject to change without notice.

For detailed technical information visit www.biomat.it