

PRODUCT DATA SHEET

FRAME FOR MICROPLATE

1. Description

Code	Format	Colour	Material	Dimension
FRG 2	12 x 8 frame with grid single well holding	White	ABS	SBS standard

2. Features of raw material

Polymer	Acrylonitrile – Butadiene - Styrene
Contains	Antistatic
Density (ISO 1183)	1050 Kg/m ³
Average Weight	17.78gr ± 0.2 gr
Packaging	Q.ty/case: 200
Other information	The moulded product is safe by the physiological aspect according to the producer information.

Data Sheet subject to change without notice.

For detailed technical information visit www.biomat.it