

PRODUCT DATA SHEET

1536 well plate GLASS BOTTOM

1. Description

| Code | Plate format | Colour | Breakable |
|-------------|-------------------------------|--------|-----------|
| MCB-GB-1536 | 1536 well plates Glass Bottom | Black | No |

| | | |
|------------|----------------------------|-------------|
| Plate | Material | Polystyrene |
| | Well bottom | Flat |
| | Well colour | Clear |
| | Well shape | Square |
| Dimensions | ANSI/SLAS 1-2004 Standards | |

2. Features

| | |
|-------------------|--|
| Packaging | Qty/ pk: 5 |
| | Qty/ cs: 40 |
| Storage | Room temperature |
| Other information | No glue inclusion into wells Glass-bottom plates are available for laser scanning and confocal imaging analysis |

Data Sheet subject to change without notice.

For detailed technical information visit www.biomat.it