

## PRODUCT DATA SHEET

### 384 well plate GLASS BOTTOM

#### 1. Description

Code	Plate format	Colour	Breakable
MCB-GB-384	384 well plates Glass Bottom	Black	No

Plate	Material	Polystyrene
	Well bottom	Flat
	Well colour	Clear
	Well shape	Square
Dimensions	ANSI/SLAS 1-2004 Standards	

#### 2. Features

Packaging	Qty/ pk: 5
	Qty/ cs: 40
Storage	Room temperature
Other information	No glue inclusion into wells Glass-bottom plates are available for laser scanning and confocal imaging analysis

Data Sheet subject to change without notice.

For detailed technical information visit [www.biomat.it](http://www.biomat.it)