

PRODUCT DATA SHEET

96 well plate GLASS BOTTOM

1. Description

Code	Plate format	Colour	Breakable
MCB-GB-96	96 well plate Glass Bottom	Black	No

Plate	Material	Polystyrene
	Well bottom	Flat
	Well colour	Clear
	Well shape	Square
Dimensions	ANSI/SLAS 1-2004 Standards	

2. Features

Packaging	Qty/ pk: 5
	Qty/ cs: 40
Storage	Room temperature
Other information	No glue inclusion into wells Glass-bottom plates are available for laser scanning and confocal imaging analysis

Data Sheet subject to change without notice.

For detailed technical information visit www.biomat.it