

PRODUCT DATA SHEET

STREPTAVIDIN HIGH BINDING CAPACITY COATED SURFACE 12 x 75 mm tubes format

1. Description

Code	Tubes format	Tubes colour	Volume per tube
TS1275STDF-SB75/500/W	12 x 75 mm	Clear	5 mL

Tubes	Material	Polystyrene
	Bottom	Round
Dimensions	Lenght	75 mm
	Diameter	12 mm

2. Features

Surface	Streptavidin HB
Volume coating	500 µl
Blocking	All tubes are pre-blocked for low non specific binding with Blocker Well (non proteic blocking buffer)
Volume blocking	500 µl
Packaging	20 tubes are packed in a barrier bag resealable with zip with desiccants
Storage	2 – 8 °C. Once opened, place unused tubes in a resealable bag with desiccant and store at 2 – 8 °C. Tubes are shipped at room temperature
Shelf life	2 years from date of manufacture
Sterile	No
Ready to use	Yes
Other information	All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

For detailed technical information visit www.biomat.it