

PRODUCT DATA SHEET

STREPTAVIDIN COATED SURFACE

PCR strip tubes format

1. Description

Code	Tubes format	Tubes colour	Volume
PCR0STF-SA5/100	8 strip PCR tubes w/o caps	Clear	0.2 mL

Tubes	Material	Polypropylene
	Bottom	Round

2. Features

Surface	Streptavidin
Volume coating	100 µl
Blocking	All tubes are pre-blocked for low non specific binding
Volume blocking	200 µl
Biotin binding	~ 5 pmol/tube (100 µl)
Packaging	10 tubes are packed in a barrier bag resealable with zip with desiccant
Storage	2 – 8 °C. Once opened, place unused tubes in a resealable bag with desiccant and store at 2 – 8 °C. PCR strip tubes are shipped at room temperature
Shelf life	2 years from date of manufacture
Sterile	No
Ready to use	Yes
Other information	Strip tubes are compatible with standard 0.2 mL thermal cycler blocks. Separate flat caps are available in strips of 8; they are easy to open and close without tools. All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it