

PRODUCT DATA SHEET

TISSUE CULTURE SURFACE

24 well plate - 96 well plate – 384 well plate

1. Description

Code	Plate format	Volume per well	Sterile	Lids
MC23FL-24-TC/S	24 well solid plate	3.4 ml	Yes	Yes
MC23FL-24-TC	24 well solid plate	3.4 ml	No	No
MC23FL-96-TC/S	96 well solid plate	400 µl	Yes	Yes
MC23FL-96-TC	96 well solid plate	400 µl	No	No
MC23FL-384-TC/S	384 well solid plate	120 µl	Yes	Yes
MC23FL-384-TC	384 well solid plate	120 µl	No	No

Plate	Material	Polystyrene
	Well bottom	Flat
Dimensions	SBS standard	

2. Features

Surface	Tissue culture treated
Active Surface	see "Volume per well"
Packaging	Plate are packed in a transparent bag
Storage	Room Temperature
Shelf life	N/A
Ready to use	Yes

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it