

## PRODUCT DATA SHEET

### MALEIMIDE COATED SURFACE

PCR strip tubes format

#### 1. Description

| Code      | Tubes format               | Tubes colour | Volume |
|-----------|----------------------------|--------------|--------|
| PCR22-MAL | 8 strip PCR tubes w/o caps | Clear        | 0.2 mL |

  

|       |          |               |
|-------|----------|---------------|
| Tubes | Material | Polypropylene |
|       | Bottom   | Round         |

#### 2. Features

|                   |   |
|-------------------|---|
| Surface           | Maleimide   |
| Volume coating    | 100 µl  |
| Blocking          | All tubes are pre-blocked for low non specific binding with non proteic blocking buffer   |
| Volume blocking   | 200 µl  |
| Packaging         | 10 tubes are packed in a barrier bag resealable with zip with desiccant   |
| Storage           | 2 – 8 °C.<br>Once opened, place unused tubes in a resealable bag with desiccant and store at 2 – 8 °C. PCR tubes are shipped at room temperature.   |
| Shelf life        | 3 months from date of manufacture   |
| Sterile           | No  |
| Ready to use      | Yes   |
| Other information | Strip tubes are compatible with standard 0.2 mL thermal cycler blocks.<br>Separate flat caps are available in strips of 8; they are easy to open and close without tools.<br>All lots are tested for uniformity and reproducibility<br>Certificate of Quality is released for every lot |

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit [www.biomat.it](http://www.biomat.it)