

PRODUCT DATA SHEET

CARBOXYLATED SURFACE PCR strip tubes format

1. Description

Code	Tubes format	Tubes colour	Volume
PCR04-COOH	8 strip PCR tubes w/o caps	Clear	0.2 mL

Tubes	Material	Polypropylene
	Bottom	Round

2. Features

Surface	Carboxylated
Active Surface	200 µl
Packaging	10 tubes are packed in a barrier bag resealable with zip
Storage	Room Temperature. Once opened, place unused strips or plate in a transparent bag and store at room temperature.
Shelf life	2 years from date of manufacture
Sterile	No
Ready to use	Yes
Other information	Strip tubes are compatible with standard 0.2 mL thermal cycler blocks. Separate flat caps are available in strips of 8; they are easy to open and close without tools. All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it