

PRODUCT DATA SHEET

CARBOXYLATED SURFACE

PCR plate format

1. Description

Code	Plate format	Colour	Volume per tubes
PCR960104-COOH	PCR plate, 96 well, non-skirted	Clear	0.2 mL
PCR960204-COOH	PCR plate, 96 well, semi-skirted	Clear	0.2 mL
PCR960304-COOH	PCR plate, 96 well, skirted	Clear	0.2 mL

Tubes	Material	Polypropylene
	Bottom	Round

2. Features

Surface	Carboxylated
Active Surface	200 µl
Packaging	Each plate is packed in a barrier bag resealable with zip
Storage	Room Temperature.
Shelf life	2 years from date of manufacture
Sterile	No
Ready to use	Yes
Other information	PCR plates are compatible with standard 0.2 mL thermal cycler blocks. Each plate is alphanumeric coding All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it