

PRODUCT DATA SHEET

CARBOXYLATED SURFACE

PCR plate format

1. Description

Code	Plate format	Colour	Volume per tubes
PCR960104-COOH	PCR plate, 96 well, non-skirted	Clear	0.2 mL
PCR960204-COOH	PCR plate, 96 well, semi-skirted	Clear	0.2 mL
PCR960304-COOH	PCR plate, 96 well, skirted	Clear	0.2 mL

Tubes	Material	Polypropylene
	Bottom	Round

2. Features

Surface	Carboxylated
Active Surface	200 µl
Packaging	Each plate is packed in a barrier bag resealable with zip
Storage	Room Temperature.
Shelf life	2 years from date of manufacture
Sterile	No
Ready to use	Yes
Other information	<p>PCR plates are compatible with standard 0.2 mL thermal cycler blocks.</p> <p>Each plate is alphanumeric coding</p> <p>All lots are tested for uniformity and reproducibility</p> <p>Certificate of Quality is released for every lot</p>

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it