

## PRODUCT DATA SHEET

### MALEIMIDE COATED SURFACE

#### PCR plate format

#### 1. Description

Code	Plate format	Colour	Volume per tubes
PCR960122-MAL	PCR plate, 96 well, non-skirted	Clear	0.2 mL
PCR960222-MAL	PCR plate, 96 well, semi-skirted	Clear	0.2 mL
PCR960322-MAL	PCR plate, 96 well, skirted	Clear	0.2 mL

Tubes	Material	Polypropylene
	Bottom	Round

#### 2. Features

Surface	Maleimide
Volume coating	100 µl
Blocking	All plates are pre-blocked for low non specific binding with non proteic blocking buffer
Volume blocking	200 µl
Packaging	Each plate is packed in a barrier bag resealable with zip with desiccant
Storage	2 – 8 °C. Plates are shipped at room temperature.
Shelf life	3 months from date of manufacture
Sterile	No
Ready to use	Yes
Other information	PCR plates are compatible with standard 0.2 mL thermal cycler blocks. Each plate is alphanumeric coding All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit [www.biomat.it](http://www.biomat.it)