

PRODUCT DATA SHEET

MOUSE MONOCLONAL ANTIBODY ISOTYPING COATED SURFACE 96 well plate strip format

1. Description

Code	Plate format	Strips colour	Volume per well	Breakable
MG36F-MAI	8 well strip on single well holding frame	Clear	350 µl	Yes
Strip	Material	Polystyrene		
	Well bottom	Flat		
	Well shape	Round		
Frame	Colour	White		
	Material frame	ABS		
Dimensions	SBS standard			

2. Features

Surface	Antibodies specific for mouse IgG1, IgG2a, IgG2b, IgG3, IgA, IgM, kappa light chain and lambda light chain
Volume coating	100 µl
Blocking	All plates are pre-blocked for low non specific binding
Volume blocking	200 µl
Uniformity	CV% < 10% with our standard test. For further details see Certificate of Quality
Packaging	Each plate is packed in a barrier bag resealable with zip with desiccant
Storage	2 – 8 °C. Once opened, place unused strips or plates in a resealable bag with desiccant and store at 2 – 8 °C. Plates are shipped at room temperature
Shelf life	24 months from date of manufacture
Sterile	No
Ready to use	Yes
Other information	All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Also available for luminescence and fluorescence assays.

Other formats are available upon request.

For detailed technical information visit www.biomat.it