

PRODUCT DATA SHEET

96 Deep Well Plate 2.2 ml
 KingFisher™ compatible

1. Description

Code	Plate format	Colour	Volume per well
DWSP2200KF	96 Deep Well plate KingFisher™ compatible	Clear	2.2 ml

Plate	Material	Polypropylene (high quality)
	Well shape	Square
	Well bottom	Pyramid
Dimensions	127 mm x 85 mm x 45 mm (L x W x H)	

2. Features

Packaging	Qty/ pk: 5
	Qty/ cs: 50
Storage	Room temperature in a dry place. Keep out of sunlight.
Other information	DNase and RNase free Good chemical resistance and heat-resistance. For use with KingFisher™ Flex, Apex and Presto.

Data Sheet subject to change without notice.

For detailed technical information visit www.biomat.it