

PRODUCT DATA SHEET

PCR 8 STRIP TUBES

1. Description

Code	Format	Colour	Volume per tubes
PCRS8	PCR strip tubes – 8 wells w/o caps	Clear	0.2 mL

Tubes	Material	High Quality Polypropylene
-------	----------	----------------------------

2. Features

Appearance	each tube is numerically identified (from 1 to 8)
Tubes bottom	Round
Tubes shape	Conical
Temperature range	- 80 °C / + 121 °C
Packaging	Qty/ cs: 125
Storage	Room temperature in a dry place. Keep out of sunlight.
Other informations	For research use only Good chemical resistance Good heat resistance Free of DNase, RNase and endotoxin Free of DNA

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it