

PRODUCT DATA SHEET

PCR PLATE

1. Description

Code	Plate format	Cut corner	Colour	Volume per tubes
PCR960102	PCR plate, 96 well, non-skirted	H1	Clear	0.2 mL
PCR960201	PCR plate, 96 well, semi-skirted	H12	Clear	0.1 mL
PCR960201A	PCR plate, 96 well, semi-skirted ABI type	A1	Clear	0.1 mL
PCR960202	PCR plate, 96 well, semi-skirted	A12	Clear	0.2 mL
PCR960301	PCR plate, 96 well, skirted	H1	Clear	0.1 mL
PCR960302	PCR plate, 96 well, skirted	H1	Clear	0.2 mL

Tubes	Material	High Quality Polypropylene
-------	----------	----------------------------

2. Features

Appearance	alphanumeric coding for optimal sample identification
Temperature range	- 80 °C / + 121 °C
Packaging	Qty/ cs: 50
Storage	Room temperature in a dry place. Keep out of sunlight.
Other informations	For research use only Good chemical resistance Good heat resistance (can be frozen for long time) Free of DNase, RNase and endotoxin Free of DNA

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it