

PRODUCT DATA SHEET

POLY-D-LYSINE COATED SURFACE

96 well plate with lid

1. Description

Code	Plate format	Colour	Volume per well
MC24F-LYS-D/A	96 well solid plate	Clear	400 µl
MCW24F-LYS-D/A	96 well solid plate	White	400 µl
MCB24F-LYS-D/A	96 well solid plate	Black	400 µl
MCWT24F-LYS-D/A	96 well solid plate	White with transparent bottom	350 µl
MCBT24F-LYS-D/A	96 well solid plate	Black with transparent bottom	350 µl

Plate	Material	Polystyrene
	Well bottom	Flat
	Well shape	Round
Dimensions	SBS standard	

2. Features

Surface	Poly-D-lysine (MW: 70 – 150 kDa)
Approx. growth area	0.32 cm ²
Packaging	5 plates are packed in a transparent bag
Storage	2 – 8 °C. Plates are shipped at room temperature
Shelf life	2 years
Sterile	produced in aseptic conditions
Lid	Yes
Ready to use	Yes

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it