

PRODUCT DATA SHEET

POLY-D-LYSINE COATED SURFACE

96 well plate – sterile with lid

1. Description

| Code | Plate format | Colour | Volume per well |
|-----------------|---------------------|-------------------------------|-----------------|
| MC24F-LYS-D/S | 96 well solid plate | Clear | 400 µl |
| MCW24F-LYS-D/S | 96 well solid plate | White | 400 µl |
| MCB24F-LYS-D/S | 96 well solid plate | Black | 400 µl |
| MCWT24F-LYS-D/S | 96 well solid plate | White with transparent bottom | 350 µl |
| MCBT24F-LYS-D/S | 96 well solid plate | Black with transparent bottom | 350 µl |

| | | |
|------------|--------------|-------------|
| Plate | Material | Polystyrene |
| | Well bottom | Flat |
| | Well shape | Round |
| Dimensions | SBS standard | |

2. Features

| | |
|---------------------|---|
| Surface | Poly-D-lysine (MW: 70 – 150 kDa) |
| Approx. growth area | 0.32 cm ² |
| Packaging | 5 plates are packed in a transparent bag |
| Storage | 2 – 8 °C. Plates are shipped at room temperature |
| Shelf life | 2 years |
| Sterile | Yes |
| Lid | Yes |
| Ready to use | Yes |

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it