

PRODUCT DATA SHEET
POLY-L-LYSINE COATED SURFACE
 384 well plate – sterile with lid

1. Description

Code	Plate format	Colour	Volume per well
MC42F-384-LYS-L/S	384 well solid plate	Clear	120 µl
MCW42F-384-LYS-L/S	384 well solid plate	White	120 µl
MCBT42F-384-LYS-L/S	384 well solid plate	Black with transparent bottom	120 µl

Plate	Material	Polystyrene
	Well bottom	Flat
	Well shape	Round
Dimensions	SBS standard	

2. Features

Surface	Poly-L-lysine (MW: 70 – 150 kDa)
Approx. growth area	0.06 cm ²
Packaging	5 plates are packed in a transparent bag
Storage	2 – 8 °C. Plates are shipped at room temperature
Shelf life	1 year from date of manufacture
Sterile	Yes
Lid	Yes
Ready to use	Yes

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it