

## PRODUCT DATA SHEET

### TISSUE CULTURE SURFACE

24 well plate - 96 well plate – 384 well plate

#### 1. Description

Code	Plate format	Plate colour	Volume per well	Approx. growth area	Sterile	Lids
MCWT23FL-24-TC/S	24 well plate	White with transparent bottom	3.4 ml	1.9 cm <sup>2</sup>	Yes	Yes
MCBT23FL-24-TC/S	24 well plate	Black with transparent bottom	3.4 ml	1,9 cm <sup>2</sup>	Yes	Yes
MC23FL-96-TC/S	96 well plate	Clear	400 µl	0.32 cm <sup>2</sup>	Yes	Yes
MCW23FL-96-TC/S	96 well plate	White	400 µl	0.32 cm <sup>2</sup>	Yes	Yes
MCB23FL-96-TC/S	96 well plate	Black	400 µl	0.32 cm <sup>2</sup>	Yes	Yes
MCWT23FL-96-TC/S	96 well plate	White with transparent bottom	350 µl	0.32 cm <sup>2</sup>	Yes	Yes
MCBT23FL-96-TC/S	96 well plate	Black with transparent bottom	350 µl	0.32 cm <sup>2</sup>	Yes	Yes
MC23FL-384-TC/S	384 well plate	Clear	120 µl	0.06 cm <sup>2</sup>	Yes	Yes
MCW23F-384-TC/S	384 well plate	White	120 µl	0.06 cm <sup>2</sup>	Yes	Yes
MCB23F-384-TC/S	384 well plate	Black	120 µl	0.06 cm <sup>2</sup>	Yes	Yes
MCWT23F-384-TC/S	384 well plate	White with transparent bottom	120 µl	0.06 cm <sup>2</sup>	Yes	Yes
MCBT23F-384-TC/S	384 well plate	Black with transparent bottom	120 µl	0.06 cm <sup>2</sup>	Yes	Yes

Plate	Material	Polystyrene
	Well bottom	Flat
Dimensions	SBS standard	

#### 2. Features

Surface	Tissue culture treated
Active Surface	see table
Packaging	bag of 4 plates / case of 48 plates (24 well format) bag of 5 plates / case of 50 plates (other formats)
Storage	Room Temperature
Shelf life	N/A
Ready to use	Yes

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit [www.biomat.it](http://www.biomat.it)