

PRODUCT DATA SHEET

AMMINATED SURFACE

96 well plates

1. Description

Code	Plate format	Strips colour	Volume per well	Breakable
MG02F-AM1	8 well strip on single well holding frame	Clear	350 µl	Yes
MGW02F-AM1	8 well strip on single well holding frame	White	350 µl	Yes
MGB02F2-AM1	8 well strip on single well holding frame	Black	360 µl	Yes
MT02F4-AM1	8 well strip on 12 x 8 frame	Clear	360 µl	No
MTW02F4-AM1	8 well strip on 12 x 8 frame	White	360 µl	No
MTB02F4-AM1	8 well strip on 12 x 8 frame	Black	360 µl	No
MC02F-AM1	96 well solid plate	Clear	400 µl	No
MCW02F-AM1	96 well solid plate	White	400 µl	No
MCB02F-AM1	96 well solid plate	Black	400 µl	No

Strip/plates	Material	Polystyrene
	Well bottom	Flat
	Well shape	Round
Frame	Colour	White
	Material frame MT	Polypropylene + glass fiber
	Material frame MG	ABS
Dimensions	SBS standard	

2. Features

Surface	Primary Amino Groups
Active Surface	200 µl
Blocking	N/A
Volume blocking	N/A
Packaging	Plates are packed in a transparent bag
Storage	Room Temperature. Once opened, place unused strips or plates in a transparent bag and store at room temperature.
Shelf life	2 years from date of manufacture
Sterile	No
Ready to use	Yes
Other information	All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit www.biomat.it