

PRODUCT DATA SHEET

NO BINDING SURFACE

96 well plate strip format

1. Description

Code	Plate format	Strips colour	Volume per well	Breakable
MG03F- NB	8 well strip on single well holding frame	Clear	350 µl	Yes
MGW03F- NB	8 well strip on single well holding frame	White	350 µl	Yes
MGB03F2- NB	8 well strip on single well holding frame	Black	360 µl	Yes
MT03F4-NB	8 well strip on 12 x 8 frame	Clear	360 µl	No
MTW03F4- NB	8 well strip on 12 x 8 frame	White	360 µl	No
MTB03F4- NB	8 well strip on 12 x 8 frame	Black	360 µl	No

Strip	Material	Polystyrene
	Well bottom	Flat
	Well shape	Round
Frame	Colour	White
	Material frame MT	Polypropylene + glass fiber
	Material frame MG	ABS
Dimensions	SBS standard	

2. Features

Surface	No Binding
Active Surface	300 µl
Blocking	N/A
Volume blocking	N/A
IgG binding	None
Uniformity	CV% < 10% with our standard test. For further details see Certificate of Quality
Packaging	Plates are packed in a transparent bag
Storage	Room Temperature. Once opened, place unused strips or plates in a transparent bag and store at room temperature.
Shelf life	3 years from date of manufacture
Sterile	No
Ready to use	Yes
Other information	All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot

Data Sheet subject to change without notice.

Other formats are available upon request.

For detailed technical information visit <https://www.biomat.it/>